

Nuclear Imaging I (SPECT)

RT4220 - LECTURE #12


KEDREE PROFFITT

WSU

10/27/2025

Basics of Nuclear Medicine

- Patient presents with a disease → suspected disease can be imaged with a particular radioactive tracer → radioactive tracer is injected → an image is formed of the disease by a machine
- The underlined section is the only difference between **gamma cameras**, single photon emission computed tomography (**SPECT**), positron emission tomography (**PET**), and now Compton cameras
- All nuclear medicine detectors are either scintillator or semiconductor based
- Nuclear imaging **permits functional imaging**, exceptionally important!



Modes of Nuclear Imaging I


- Gamma Camera
 - Planar images via projections
 - Typically scintillators and PMTs
 - Physical collimation
 - Uses low-energy photon emitters
- SPECT
 - Tomographic imaging (3D) via step-and-shoot planar images
 - Typically scintillators and PMTs
 - Physical collimation
 - Uses low-energy photon emitters

Modes of Nuclear Imaging II

- PET
 - Intrinsically tomographic data (3D)
 - Typically scintillators, but moving towards semiconductors
 - Electronic collimation
 - Uses positron emitters (high-energy photons)
 - Generally better spatial resolution and sensitivity than SPECT
 - Functional imaging

Gamma Camera Basics


1. Injection
2. Gamma photons generated by beta- emitter or Tc-99m
3. Photons leave patient and pass through collimator
4. Remaining photons are scintillated (each keV produces roughly 38 low energy photons in the visible band for NaI(Tl))
5. PMT converts scintillated photons into electronic signal
6. A pulse height analyzer rejects suspected non-valid events
7. A position logic circuit determines the origin location of each signal, now it is counted!
8. Many counts are summed together to form a planar image









TOMOGRAPHIC IMAGING SPECT (SPECT/CT)

- Acquisition of multiple planar images in a **step-and-shoot** fashion from **various angles** around the patient allows 3D image recon → **64 projections 10-20s each for 180 arc (128 projections for 360 arc)**




Step Three: Collimation


- Parallel hole is most common in nuc. med.
- Pinhole is used for thyroid uptake analysis
- There are other forms of collimation that are gaining prominence, but are too complicated for us




Step Four: Possible Photon Fates




Step Six: Discrimination of Energy



Step Seven: Event Localization



Step Seven Continued



Common Procedures

- Ask the nuc. med. people what they see most (Almost all are Tc-99m based)
- Bone Scans
 - Cancer
 - Occult Fractures
 - Osteomyelitis/Arthritis
- Cardiac Scans
- Metastatic Cancer
- Prostate Cancer
- Neurological Scans

INTRINSIC / EXTRINSIC RESOLUTION


Intrinsic



Detector

Crystal

Bar phantom



Point Source

Extrinsic



Flood Source


Collimator

PERFORMANCE SPECIFICATION

- **Intrinsic** spatial resolution (3.5 – 4.5 mm)
- System (**extrinsic**) spatial resolution (8 - 12 mm)
- **Energy** resolution (9 – 10%)
- **Intrinsic sensitivity** 200-300 cpm/uCi (~ **0.7%**)
- **Clinical sensitivity** (number of photons generated in the patient relative to the number of photons measured) is only ~ 1 in 10,000 (**0.01%**)
- This is due to absorption/**Compton scatter** in patient and absorption by the **collimator** (to create image)


GAMMA CAMERA QUALITY CONTROL

- Daily QC > Flood phantom to test uniformity
(Differential / Integral Uniformity < 5%)
Peaking of energy (remove background)
- Weekly QC > Bar phantom to monitor spatial resolution
- Monthly QC > Jaszczak phantom to test spatial resolution
FWHM ~ 1.7 x smallest pattern
- Quarterly QC > Sensitivity test
- Annual QC > Rotational uniformity
Dead-time test




bar
phantom

Jaszczak
phantom




SPECT CORRECTIONS

- Regional spectra shifts > **Energy Correction**
- Miss-positioning of events > **Spatial Linearity Correction**
- Regional count differences > **Uniformity Correction**
- Count rate losses > **Pulse Pile-up Correction**



Energy correction




Linearity correction

NEW DEVELOPMENTS IN NUCLEAR MEDICINE

- Novel semiconductor detectors (CZT)
- SPECT Ring design using Focused Detectors
- Use of SiPM instead of PMT for scintillation crystals
- Photon counting technology



CADMIUM ZINC TELLURIDE (CZT) SEMICONDUCTOR DETECTOR

- **Direct** conversion of gamma rays into electrical charge
- **Room temperature** semiconductor (E gap ~ 1.6 eV) with high resistance resulting in very **low leakage** current
- High **sensitivity** due to high density (5.8 g/cm^3 , NaI - 3.7 g/cm^3)
- Excellent E resolution ($\sim 5\%$, NaI $>10\%$)



SPECT/CT


- Conventional SPECT uses a **Two-head** detector system in **step-and-shoot mode**
- **12 digital CZT detectors** that automatically move within mm of the patient's body




SPECT RING DESIGN

- Ring design allows **faster scanning** due to increased detector number and **decreased distance** between detector and patient

Conventional design



Ring design



SPECT RING DESIGN ADVANTAGES / DRAWBACKS

- **Advantages**

- Detectors closer to patient – improved resolution
- No “wasted” detector space
- Only 2-6 steps required – faster scanning
- CZT semiconductor detector – compact detector

- **Drawbacks**

- Non-exchangeable collimator – limited energy range
- No true sinogram
- Expensive

LEHR collimator:
40 – 250 keV

FROM PMT TO PHOTON COUNTING USING SiPM

- **APD** → Avalanche Photodiode provides amplification **prior** to breakdown (**analog signal**)
- **SiPM** → array of small APDs, each of which is in **Geiger mode** allows photon counting (**digital signal**)




	PMT	APD	Analog SiPM	Digital photon counting Digital SiPMs
TOF capability*	● analog, time-of-flight	● N/A	● analog, time-of-flight	● digital, time-of-flight
Operational stability	● medium	● low	● medium	● high
Signal amplification	● 10^6	● 10^{2-3}	● 10^6	● not needed
Level of Integration	● low	● medium	● medium	● high
Signal readout	● analog	● analog	● analog	● digital

AVALANCHE PHOTO DIODE (APD)

Voltage **below** breakdown = **Avalanche** mode

Voltage **above** breakdown = **Geiger** mode



MODE OF OPERATION DEPENDS ON BIAS VOLTAGE

- Analog to gas detectors (but “**Avalanche**” is used differently)




Photo Diode (PD)

Only primary electrons

Avalanche PD (APD)

Linear amplification
 $V(bias) < V(breakdown)$

Single Photon APD


Continuous current
 $V(bias) > V(breakdown)$

ANALOG SILICON PM (SiPM)

2D array of APDs, biased above breakdown voltage (Geiger mode)

Silicon Photo-Multiplier (SiPM)

- Compact Device
- Operating voltage (30-120V)
- Resolution - Single photon detection
- Response time – ~100 ps
- High gain - 10^6
- High Quantum Efficiency – 90%
- High Photon Detection Efficiency – 60%
- Immunity to Magnetic Field



Each cell 25 x 25 μm^2

