



Nuclear Imaging II (PET)

RT4220 - LECTURE #13


KEDREE PROFFITT

WSU

10/30/2025

Basics of PET Imaging

- Positron Emission Tomography → Positron → 2x 511 keV photons oppositely opposed → *near*-simultaneous detection of each photon tells us where it came from → we can **quantify** how much tracer was deposited in a given location
- Extremely short half-life of positron emitters:
 - C-11: 20 mins
 - O-15: 2 mins
 - F-18: 110 mins
- Short half life **requires** on-site production; we have a generator!
- Most common PET radiopharmaceutical is FluoroDeoxyGlucose (F-18 != FDG)




LSO Detector
ring

Vision 600 PET/CT
TOF = 213 ps


Rotating CT
anode

Continuous
bed motion



FDG Scans


- Types of Scans:
 - Tumor imaging and staging
 - Detection of sub-clinical involvement
 - Cardiac tissue viability studies
 - Brain Imaging
- Function:
 - Localizes like glucose in the body but gets stuck




CT

PET

Fusion




CT Image



Poorly Defined Tumor Margins

Fused CT-PET Image



FDG Avid Tumor

Common PET Usage

- Primarily to detect cancer and monitor its progression/metastasis status
- PET/CT has become the primary **imaging** tool for cancer staging
 - There are other great non-imaging tools!
- Apparently, some are using PET to diagnose Alzheimer's?



Common SPECT Usage

- SPECT primarily for tracking the progression of heart disease (blocked coronary arteries)
- There are also radiotracers for detecting bone, gall bladder, and intestinal disorders
- Apparently, some people have used it to help diagnose Parkinson's?


PET ALLOWS AUTORADIOGRAPHY **IN VIVO**

- Positron Emission Tomography allows **radiography in vivo** - **absolute quantification** of **radioactivity concentration** in different parts of the human body (due to **exact** attenuation correction)
- Use of **electronic collimation** instead of lead collimation - high sensitivity and efficiency (**picomolar** concentration range)
- Use of **short-lived** biologically active radio-pharmaceutical: C-11flumazenil, N-13 ammonia, O-15 water, F-18 deoxy-glucose

PET DATA ACQUISITION AND ANALYSIS



FDG workflow



Irradiation:
1- 2 h

Starting activity:
3-12 Ci



Assembly and
radiosynthesis:
~ 2 h

FPV:
2 - 10 Ci




Quality Control
& drug release

~ 30 min

Overall time: 3 - 5 h
Expiration: 12 h

ENERGY DISTRIBUTION OF EMITTED POSITRONS







Energy spectrum of a
positron emitted from ^{18}F



The positron travels a distance before it annihilates **dependent** on its energy

SPATIAL DISTRIBUTION OF ANNIHILATION EVENTS




Pathlength of positron is a random walk

Exponential distribution


Monte-Carlo simulation
(100,000 histories)

COINCIDENCE DETECTION

- Coincidence events occur during time window τ
 τ typically 5 ns



TIME-OF-FLIGHT (TOF) DATA ACQUISITION



Conventional:

Detected event projected to **all** voxels between detector pairs

Statistical noise from activity in one voxel adds noise to **all** voxels

⇒ Large Noise Amplification


Time-of-Flight:

Detected event projected **only** to voxels consistent w/ measured time

Statistical noise from one voxel adds noise to **only a few nearby** voxels


⇒ Small Noise Amplification

REFORMATTING THE DATA INTO SINOGRAMS




RECONSTRUCTION FROM PROJECTIONS

- Ring of detectors acquires **360 degrees** projections of **511 keV photons emitted** from patient




SINOGRAM OF A POINT SOURCE


- Path of a point source follows a **sine** curve




TYPES OF COINCIDENCE EVENTS



True



Scatter



Randoms

—————→ Gamma ray
----- Assigned LOR

Other Nuclear Imaging



ALL THE THINGS I NEVER SAW

Nuclear Medicine Therapy

- I-131 Radionuclide Therapy
 - Iodine collects in the thyroid gland
 - Small doses can induce cancer in the thyroid, large doses can kill thyroid cancers
 - Accidental exposure can be flushed out of system with non-radioactive iodine (I-127) such as potassium iodide (KI)

The thyroid controls many body processes including heart rate, blood pressure, body temperature and childhood growth and development

I-131 Therapy

- I-131 (administered orally) is intended for:
 - Treating thyroid diseases (ie. hyperthyroidism: 4-10mCi)
 - Treating thyroid carcinomas (100-150mCi)
 - Post thyroidectomy or thyroid lobectomy (75-150mCi)
 - Treatment of distant metastases (150->200mCi)
 - Diagnostic studies (I-131: 4-10μCi or I-123: 100-200 μCi)
- Therapeutic patients:
 - Must avoid contact with pregnant women, infants and young children for a few days post treatment
 - Should not sleep in the same bed with another person for a few days post treatment

Hospitals may quarantine therapeutic patients

I-131 Therapy


*Hospital room setup for a
therapeutic I-131 patient*



I-131 Therapy

- Beta minus decay
 - 8.01 day half life
 - 90% of Decays have:
 - Max beta particle energy = 606.3keV
 - Average energy of 190keV ($\sim 1/3$ MaxE)
 - Beta energy absorbed within 0.6 to 2.0mm of tissue thickness
- ~637 keV & 365 keV gamma emitted
 - Can be used for imaging or uptake measurement
 - I-131 is never used solely for diagnostic purposes

I-131 Decay Chart




I-123 Diagnostic Studies

- I-123 has a short half-life (13.3hr)
 - I-131 half-life was 8.01 days
 - 97% Beta Minus, 3% EC
 - Emits 159keV gammas
 - Allows for much lower dose to patient
- Patient studies are performed 6-24hrs after ingestion
 - Allows time for accumulation in the thyroid
- Gamma camera's NaI detector is efficient due to the 159keV photons that emitted



I-123 Diagnostic Studies

- Thyroid uptake study:
 - Determine thyroid functionality by counting how much radioactivity is taken into the thyroid
 - Normal uptake = 15-25%
 - Hypo: <10%
 - Hyper: >30%
 - Performed using a gamma camera with a **pinhole collimator** or a scintillation detector called a "thyroid probe"



I-123 Diagnostic Studies

- Imaging with Gamma Camera
 - Used to find abnormalities in absorption uptake of iodine in the thyroid



- Thyroid Uptake Probe
 - Scintillation counting system
 - Used to find uptake of iodine in the thyroid
 - Images are not created with the uptake probe



Iodine Thyroid Uptake Study

- Radioiodine thyroid uptake:


$$\% \text{ *Uptake* } = \frac{\text{Neck counts (cpm)} - \text{Thigh Counts(cpm)}}{\text{Administed counts (cpm)} - \text{Background (cpm)}} \times 100\%$$

- The “administered counts” are found by measuring the ingested dose in a neck phantom and correcting for decay
- The “thigh counts” subtract out background absorption in the body




Theranostics

- Thera(-peutic / Diag-)nostics
- Links radionuclides to specific receptors
 - SSTR2 – Somatostatin Receptor
 - Expressed in NETs
 - NeuroEndocrine Tumors
- Ga-68 is a positron emitter (PET)
- Lu-177 DOTATATE = “Lutathera”



Current Theranostic Treatments of Interest

- Xofigo: Ra-223 Dichloride
 - Bone Metastases following Prostate Cancer
- Unnamed: Lu-177 PSMA
 - Prostate-Specific Metabolic Antigen
- Lutathera
 - GEP (Gastroentero-pancreatic)-NETs
- SIR-Spheres & TheraSphere: Y-90 MicroSpheres
 - HCC (Hepatocellular Carcinoma, a liver cancer)
- Ac-224 PSMA
 - Alpha emitter




Trends in Theranostics

- What kind of tumors?
- GEP-NET's (small percentage)
- Prostate
 - You'll probably have to work with these someday

Typical workflow in Theranostics

- Radiation Safety Protocol Development
 - Transportation of drug from hot lab to clinic to storage
 - Clinic safety
 - Policies for the room
 - Training of staff (physician through housekeeping)
 - Patient education paperwork
- Package Acceptance, Verification, and Storage
 - Wipe test, check U/mCi, label for tx
- Preparation of Drug for Delivery
- In-Clinic Drug Delivery Assistance
- Cleanup of Room, Radioactive Waste
- Storage and Disposal of Radioactive Waste



Dosimetry in Theranostics

- Xofigo: 1.49 uCi/kg of patient weight sent, delivered in syringe
- Lu-177 PSMA: Measure out 200mCi in a 10cc syringe.
- Lutathera: Standard dose of 200mCi sent in vial, delivered to all patients via gravity method
 - Or, if you're a maverick, via 60cc syringe and a chemotherapy pump

