



# Detection of Radiation I

RT4220 – Lecture #10

Kedree Proffitt

WSU

10/21/2025

# Purpose of Radiation Detectors



We cannot sense radiation with our bodies, full-stop



We need to protect people and quantify ways to reduce dose



We need to make sure we are giving a safe dose during imaging and therapy

# Principle of Radiation Detection

Relevant modes of radiation interaction:

- Excitation (Scintillators)
- **Ionization (Gas-Filled)**

It all comes back to ionization → DNA damage

So we should measure ionizations! But we cannot stick a probe inside our bodies

We must either extrapolate or translate ionizations in air to human tissue

# Types of Radiation Detectors


- There are so many:
  - **Gaseous Detectors**
    - Dose Calibrators, Ion Chambers, GM Counters, Proportional Counters...
  - **Semi-conductor Detectors**
    - Ge, Si, Diodes
  - **Film**
    - Radiographic
    - Radiochromic
  - **Scintillation Detectors**
    - NaI, CsI, BGO, TLDs, OSLDs
  - And many more on the way



# Gaseous Detectors

Firefighters? Ok sure

# Gaseous Detector Operation



# Ionization of Gas

- The **average** energy required to ionize an atom of 'normal, dry, atmospheric gas' is:

$$\left(\frac{\bar{W}}{e}\right) = 33.97 \frac{\text{eV}}{\text{ion pair}}$$

- If there is a voltage applied to the gas, then the ion pairs will be ripped apart
- Small Voltage → ions have enough time to recombine with other ions
- High Voltage → ions can gain enough energy to ionize the air by themselves!
- Some detectors are pressurized, others are open
  - Open-air chambers must correct for  $pV = nRT$  stuff




# The Gas Amplification Curve

Six regions of operation depending on desired measurement:

- I. Not useful
- II. Ionization Chambers
- III. Proportional Counters
- IV. Not useful
- V. GM Counters
- VI. Not Useful

Gaseous detectors count either **ionizations events** (cps), or **energy deposited** (dose).





---

Low voltage ( $\sim 0\text{-}100\text{ V}$ )

## **I. Recombination Region**

---

Ion pairs are **not pulled apart and collected fast enough**  $\rightarrow$  recombination

---



**Lower than expected and unreliable** ion pair counts  $\rightarrow$  **not used**

## II. Ionization Region

- Voltage is high enough for us to collect ~all of the ion pairs (~ 50-300 V)
- Used for **clinical ionization chambers** as the signal is directly proportional to ionizations (flat/straight line)
- Most relevant ion chambers are: **Farmer** (Cylindrical, Photons) and **Markus** (Plane-Parallel, Electrons)
- They tell us accurately **how much dose** a point and material is receiving



# Ionization Chambers



# III. Proportional Region

- Higher voltage than ionization region (~ 300-800 V)
- Charge collected is proportional to **energy and type of radiation**
- Gives rise to the proportional counter, it can tell us *what* type of radiation is out there



## IV. Limited Proportionality Region

High voltage (~ 800-1000 V)

```
graph TD; A[High voltage (~ 800-1000 V)] --> B["Response to collected energy diminishes while response to applied voltage increases."]; B --> C[Not used];
```

“Response to collected energy diminishes while response to applied voltage increases.”


Not used

# V. Geiger-Müller Region

- Highest useful voltage (~ 1000-1400 V)
- One ionization event is **accelerated**, creating an **avalanche** of further ionization events
- **Extremely sensitive**, but can only tell us **how many events**, nothing to do with dose
- Has to be **quenched** so that the cycle breaks and it can record a new event

## Pancake Probe Geiger Counter

Mica Window



Detachable Probe

Radiation: X-ray, Gamma, Beta, Alpha  
Operating Voltage: 900V  
Accuracy: +/-20% [ONCOLOGYMEDICALPHYSICS.COM](http://ONCOLOGYMEDICALPHYSICS.COM)

## VI. Continuous Discharge



The voltage is so high that the potential arcs through the gas, breaking stuff



**We do not want this to happen**




# Ionization Chambers

The image is a composite. The background is a blurred photograph of a person wearing a lab coat and safety glasses, looking down. In the foreground, there is a close-up of a glass pipette tip on the left and a multi-well plate on the right. The plate has many small wells, some of which contain colored liquids (red, yellow, blue).


In more depth

# Cylindrical Chambers (Farmer)

- Sometimes called a thimble chamber
- Farmer is the most common design
- 150 (QA) and (typical operation) operation modes
- Unsealed volume **300 V (correct for air pressure and temperature)**
- Standard chamber for **absolute dosimetry** of a **LINAC's output measurements**
- **Fully guarded chamber**, uses a **triaxial cable** (like coaxial, but cooler! And expensive, do not walk on them please)
- Sensitive volume of 0.6 cc
- Central electrode is made of graphite or aluminum



# Plane-Parallel Chambers (Markus)



- Markus is the most common design
- Have a **flat-top**, direction dependent
- Much **smaller sensitive volume** (0.055 cc)
- Extremely thin entrance window is more suitable for **electron measurements**

# Ionization Chamber Survey Meter

- Used to monitor radiation levels for radiation protection purposes
- Gives more information than a GM survey meter (it gives dose!)
- Battery operated and portable



# Dose Calibrator

- NOT THE SAME THING AS A WELL COUNTER, DO NOT HEAR THEIR LIES
- Used in **nuclear medicine departments** to accurately **assay patient doses**
- **Sealed** chamber eliminates  $pV=nRT$  faff
- Converts exposure to **activity** automatically
- It is an IONIZATION CHAMBER in the SHAPE OF A WELL







# Proportional Counters

In more depth!



# Proportional Counter Design

- CHERRY wants you to know that proportional counters are NOT just ion chambers at a high voltage
- Purpose of the high voltage is to allow **some acceleration ionizations**
  - Different from GM, which would continue forever without quenching! The amplification factor is lower, self-quenching
- Construction of the prop. Counter facilitates and encourages amplification
  - Uniformity
  - Filled with noble gas which allows electrons to move easier (Argon/Xenon)

# Proportional Counter Use

Can both detect **and** count individual radiation events

Size of current pulse depends on energy deposited → energy sensitive

Can distinguish between different energy events, allowing for:

- Analyzing the X-ray spectrum of an X-ray tube
- Investigating a mixed-radiation environment
- Research applications involving alpha particles



# Geiger- Muller Counters

In more depth!!!!



# GM Counter Design

- The **sealed** chamber is typically filled with argon and a quench gas like chlorine
- The quench gas soaks up excess electrons, otherwise, it would become **paralyzed**
- GM counters can only detect **one event at a time**
- If events overlap, they might only be counted as one event
- **Dead time** is how long two events need between each other to be distinguishable

## Standard Probe Geiger Counter

Detachable Probe



## Modern Pocket Geiger Counter



Multiple  
Internal  
GM Tubes  
(Mica Window)

Software  
Radiation  
Source  
Correction

Data Storage  
and Plotting

# Paralyzable Systems

- Counting systems have two categories:
  - Non-Paralyzable: If an event occurs during the deadtime of a preceding event the second event is **simply ignored**
  - Paralyzable: an event introduces a deadtime even without counting
- GM Counters **are** paralyzable!
- Classic Chernobyl scene

